

# Bottles, Plastic

## Plastic, Wide Mouth continued

### Bottles, HDPE

Wide mouth with leakproof closures in HDPE having LDPE foam inlay wads. Designed for the storage of powders. Square shape. Supplied singly.

	Ref. BGE-	Cap. ml
<b>BR800-23</b>	—	60
<b>BR800-42</b>	—	250

### Bottles, Nalgene HDPE

Wide mouth with leakproof polypropylene cap. Withstand freezer use down to -100°C.

	Ref. 2114-	Cap. ml
<b>BR807-42</b>	0008	250
<b>BR807-50</b>	0016	500
<b>BR807-80</b>	0032	1000

### Bottles, Nalgene PPCO

Polypropylene co-polymer. Wide mouth with leakproof polypropylene screw cap. Can be autoclaved at 121°C.

	Ref. 2110-	Cap. ml
<b>BR809-23</b>	0002	60
<b>BR809-34</b>	0006	175
<b>BR809-42</b>	0008	250
<b>BR809-50</b>	0016	500
<b>BR809-80</b>	0032	1000

### Bottles, HDPE

Wide mouth, square shaped, in natural or white high density polyethylene, with blue polypropylene tamper-evident cap with conical internal seal. UN-approved for the transport of solids (except 50ml sizes). Supplied singly.

	Cap. ml
<b>Natural</b>	
<b>BR815-15</b>	50*
<b>BR815-25</b>	100
<b>BR815-37</b>	250
<b>BR815-45</b>	500
<b>BR815-75</b>	100
<b>BR815-81</b>	1500
<b>BR815-94</b>	2500
<b>White</b>	
<b>BR817-15</b>	50*
<b>BR817-25</b>	100
<b>BR817-37</b>	250
<b>BR817-45</b>	500
<b>BR817-75</b>	1000
<b>BR817-81</b>	1500
<b>BR817-94</b>	2500

\*Not UN-approved

### Bottles, PVC

Transparent or amber PVC, extra-wide neck with leakproof black PP cap with PE-foam liner. Ideal for viewing contents, and easy to fill and clean. The amber versions are suitable for storing light-sensitive products. Supplied singly.

	Cap. ml	Neck i.d. mm
<b>Transparent</b>		
<b>BR865-20</b>	50	23
<b>BR865-30</b>	100	31
<b>BR865-45</b>	300	41
<b>BR865-50</b>	500	54
<b>BR865-80</b>	1000	69
<b>BR865-92</b>	2000	69
<b>Amber</b>		
<b>BR867-20</b>	50	23
<b>BR867-30</b>	100	31
<b>BR867-45</b>	300	41
<b>BR867-50</b>	500	54
<b>BR867-80</b>	1000	69
<b>BR867-92</b>	2000	69

#### Important Warning Sterilising Bottles

Some bottles where indicated can be sterilised by autoclaving. However BEFORE AUTOCLAVING (or sterilising by gas or dry heat) the cap or closure must either be removed or just set on top of the container **WITHOUT ENGAGING THE THREAD.**



BR800



BR809



BR865

BR867



BR807



BR815, BR817 are similar