

Bottles, Plastic & Culture/Specimen

Plastic, Large Capacity, Square, Octagonal and Rectangular

Bottle, HDPE, square

Natural HDPE. Ideal for storage and transport of solid materials. Space saving square design with moulded-in handgrip for safe convenient carrying. Wide neck for easy filling and supplied with blue PP cap with internal conical seal and tamper-evident ring. Overall 135 x 135 x 195mm high, neck i.d. 69mm. Supplied singly.

BR931-05 2300ml

Bottles, Azlon polypropylene, octagonal

Heavy duty rigid construction with offset neck, leakproof screwcap, approximate graduations and integral carrying handle. Can be autoclaved.

	Ref.	Capacity litres
BR932-30	BNP05B	5
BR932-40	BNP10B	10
BR932-50	BNP20B	20

Jerrycans, Azlon rigid polyethylene, rectangular

Robust, lightweight with screw cap and integral handle.

BR935-30 5 litre
BR935-50 25 litre

Jerrycans, rigid polyethylene, UN certified

Robust, stackable, food grade polyethylene jerrycans with tamper-evident screw cap.

	Capacity litres	Cap type
BR937-30	5	DIN51
BR937-40	10	DIN51
BR937-50	20	DIN61
BR937-55	30	DIN61

Aspirator stopcocks to fit BR937 jerrycans - see SV942.

Glass, Culture and Specimen

Bottles, culture

Clear glass. All caps where specified have a rubber liner.

		Capacity	Pack qty.
With aluminium screw cap, plain			
BS150-10	Bijou	7ml (¼oz)	288
BS150-25	Universal	28ml (1oz)	144
Without cap			
BS152-10	Bijou	7ml (¼oz)	255
BS152-25	Universal	28ml (1oz)	169

Caps for BS150/152

Aluminium with rubber liner. Supplied in packs of 100.

BR247-06 For 7ml and 14ml Bijou size R3/20
BR247-15 For 28ml Universal size R3/28

Caps, Sterilin, polypropylene

Wadless, non-sterile. Supplied in pack of 2850.

BS165-10 For 28ml Universal size R328

Culture vials, Universal, tubular glass

Tubular formed, Type III soda glass vials, capacity 30ml. With R3/24 neck and fitted white PP caps. Supplied in pack of 500.

BS170-25 Tubular vials, 30ml Universal size

Important Warning Sterilising Bottles

Some bottles where indicated can be sterilised by autoclaving. However **BEFORE AUTOCLAVING** (or sterilising by gas or dry heat) the cap or closure must either be removed or just set on top of the container **WITHOUT ENGAGING THE THREAD.**



BR931-05



BR935



BS150



BS170



BR932



BR937



BS152