

# Cork Products, Cotton Wool

## Cork

### Corks, tapered

Dimension given is the small end diameter. Pack quantities are approximate.

	Dia. mm	Pack qty.
<b>CW300-06</b>	6	250
<b>CW300-08</b>	8	250
<b>CW300-10</b>	10	225
<b>CW300-12</b>	12	150
<b>CW300-14</b>	14	150
<b>CW300-16</b>	16	110
<b>CW300-18</b>	18	80
<b>CW300-19</b>	19	60
<b>CW300-22</b>	22	50
<b>CW300-25</b>	25	30
<b>CW300-29</b>	29	30
<b>CW300-32</b>	32	25
<b>CW300-35</b>	35	20
<b>CW300-38</b>	38	18
<b>CW300-42</b>	42	13
<b>CW300-45</b>	45	12
<b>CW300-48</b>	48	10
<b>CW300-50</b>	50	10
<b>CW300-70</b>	Assorted	150

### Cork Rings, compressed

For supporting round bodied flasks. Depth as indicated. Supplied singly.

	o.d. x i.d. mm	For flasks ml
<b>35mm deep</b>		
<b>CW308-08</b>	80 x 30	25 to 150
<b>CW308-12</b>	110 x 60	250, 500
<b>CW308-16</b>	140 x 90	1000
<b>CW308-20</b>	170 x 120	2000, 3000
<b>CW308-24</b>	210 x 150	4000, 5000
<b>45mm deep</b>		
<b>CW309-30</b>	240 x 200	4000 to 10000

### Cork Borer Sets

Nickel plated brass. Each set complete with push rod for cleaning.

	Sets of	i.d. mm
<b>CW320-06</b>	6	4.1 to 10.3
<b>CW320-12</b>	12	4.1 to 17.8
<b>CW320-18</b>	18	4.1 to 25.3

### Cork Borer Sharpener

Suitable for all borers CW320 series.

**CW325-10** Sharpener

## Cotton Wool

For cleaning and plugging.

### Non Absorbent

In rolls of 500g.

**CW500-10** White

### Absorbent

Hospital Quality in rolls of 500g.

**CW505-10** White

Alternative Cleaning and Sealing Materials – see *Cleaning and Sealing sections.*

### Rubber Stoppers

To BS2775. In Natural, Neoprene and Silicone Rubbers. Each stopper has a moulded-in number, corresponding to the measurement of the bottom diameter in millimetres. Sizes from 3 to 105mm. 1- and 2-hole versions available in natural rubber.

– see *SW600/SW684 series in the Stopper section.*

### Stoppers, Steristoppers

Preformed cellulose plugs for sealing culture tubes and flasks. They are permeable to air and can be hot air sterilised up to 200°C. All stoppers are tapered and the cellulose material will depress slightly to ensure a tight fit into the tube or neck. Length varies from approximately 30 to 60mm according to size.

– see *SW536 series in the Stopper section.*



**CW300**



**CW308/CW309**



**CW325**



**CW320**



**CW500, CW505 is similar**