Filtration

Membrane Filter Holders and Filtration Systems

Whatman'

Vacuum Filtration Systems

Stainless steel with metal, standard or rapid-fix clamp as indicated. Maximum working temperature 200°C. Autoclavable and can be dry-heat sterilised at up to 180°C. With stainless steel screen membrane support for 47/50mm diameter membranes. Can be mounted on an accessory filter flask using an appropriate size stopper.

	Clamp	Reservoir, ml	Pre-filter size, mm	
FC304-09	Standard	500	40	
FC304-14	Rapid-fix	500	40	

Manifold Filtration Systems

For use with vacuum filter holders, as indicated. Stainless steel. With 9mm i.d. tubing connection ferrule.

FC316-10 3-place manifold, empty

3-place manifold, with 3 x FC304-09 FC316-14 FC316-16 3-place manifold, with 3 x FC304-14

FC316-30

6-place manifold, empty 6-place manifold, with 6 x FC304-09 FC316-34

FC316-36 6-place manifold, with 6 x FC304-14





