Conical

All capacities are nominal.

Flasks, Azlon, Screw Capped

Wide neck. With polypropylene screw cap and approximate graduations. Will also accept accessory conical stoppers as indicated. Autoclavable at 121°C. Supplied singly.

| Ref. | Capacity ml | Cap size mm | Stopper accepted |
|------------|---|--|--|
| translucen | t | | |
| FDE542 | 50 | 40 | 34/35 |
| FDE544 | 100 | 40 | 34/35 |
| FDE546 | 250 | 52 | 45/40 |
| FDE548 | 500 | 52 | 45/40 |
| FDE550 | 1000 | 52 | 45/40 |
| sparent | | | |
| FDE444 | 100 | 40 | 34/35 |
| FDE446 | 250 | 52 | 45/40 |
| FDE448 | 500 | 52 | 45/40 |
| FDE450 | 1000 | 52 | 45/40 |
| | translucen FDE542 FDE544 FDE546 FDE548 FDE550 sparent FDE444 FDE446 FDE448 | ml ml translucent FDE542 50 FDE544 100 FDE546 250 FDE548 500 FDE550 1000 sparent FDE444 100 FDE444 250 FDE446 250 FDE448 500 | ml mm translucent mm FDE542 50 40 FDE544 100 40 FDE546 250 52 FDE548 500 52 FDE550 1000 52 sparent FDE444 100 40 FDE446 250 52 52 FDE448 500 52 52 |

Accessory stoppers - see SW432/SW446 series

Tissue Culture

Erlenmeyer Flasks, Corning

Polycarbonate with polypropylene screw cap which can be adjusted to provide gassing or a leak-proof seal. Transparent and graduated. Sterilised by gamma irradiation. Sealed individually in a sterile sleeve within outer packs as indicated.

| | Ref. | Capacity ml | Cap size | Pack qty. | |
|----------|--------|----------------|-------------|--------------|--|
| TP550-15 | 430421 | 125 | 26 | 50 | |
| TP550-20 | 430183 | 250 | 31 | 50 | |
| TP550-25 | 430422 | 500 | 43 | 25 | |

Culture Flasks, baffled, Nalgene

Tough polycarbonate flasks with indented baffles which increase mixing when trypsinizing and aerating in a variety of processes such as cell culture growth and fermentation. Supplied in packs as indicated.

| | Ref. | Capacity ml | Pack qty. | |
|----------|-----------|----------------|--------------|--|
| TP575-13 | 4110-0250 | 250 | 4 | |
| TP575-16 | 4110-0500 | 500 | 4 | |
| TP575-19 | 4110-1000 | 1000 | 2 | |

Further ranges of these flasks can be found in the Tissue Culture section.



FK767



TP575





© 0207 729 7998 e-mail: customer.services@thelabwarehouse.com fax orders: 020 7729 2657