

# Funnels

## Buchner and Hirsch Type

All capacities and dimensions are nominal. All funnels are supplied singly.

### Buchner Funnels

Porcelain, with perforated support plate. Dimension given is the nominal filter paper size.

	Cap. ml	Paper cm
<b>FT804-12</b>	25	4.0
<b>FT804-18</b>	85	5.5
<b>FT804-26</b>	150	7.0
<b>FT804-34</b>	265	9.0
<b>FT804-40</b>	510	11.0
<b>FT804-45</b>	700	12.5

### Buchner Funnels

Polypropylene. Robust, two-piece construction for easy cleaning. Autoclavable. Dimension given is the nominal filter paper size.

	Cap. ml	Paper cm
<b>FT815-12</b>	40	4.25
<b>FT815-18</b>	70	5.5
<b>FT815-26</b>	180	7.0
<b>FT815-34</b>	390	9.0
<b>FT815-40</b>	810	11.0
<b>FT815-48</b>	2100	15.0
<b>FT815-56</b>	6000	24.0

### Vacuum Ring Seal

A neoprene rubber disc with a tapered raised boss and central hole which accepts the stem of a Buchner or similar filter funnel. Ring disc diameter 75mm.

**FK364-10** Vacuum ring seal

Buchner Rings – see FK364-25 to -75 series

### Three-piece Funnel, Whatman

Borosilicate glass with acrylic support plate for filter paper. Can be dismantled quickly for cleaning. Dimension given is the filter paper size required.

	Cap. ml	Paper cm
<b>FT850-10</b>	16	2.5
<b>FT850-15</b>	36	4.7
<b>FT850-20</b>	115	7.0
<b>FT850-22</b>	210	7.0
<b>FT850-24</b>	400	7.0
<b>FT850-30</b>	200	9.0
<b>FT850-35</b>	530	12.5

### Hirsch Funnels

Porcelain, with perforated support plate.

	Top dia. mm	Plate dia. mm	Perforated dia. mm
<b>FT880-54</b>	50	30	25
<b>FT880-58</b>	78	47	33



**FK364-10 in use**

## Thistle

### Thistle Funnels

Glass with plain stem. Stem 6mm diameter.

	Stem length mm
<b>FT922-20</b>	200
<b>FT922-25</b>	250
<b>FT922-30</b>	300

### Thistle Funnel, Azlon

Polypropylene. With stub stem. Includes stem tubing 250mm long to fit inside stub. Top diameter 34mm.

**FT930-10** FWC220

Tissue Culture bottle top filters – see FD210 series.



**FT930**



**FT804**




**FT850**



**FT815**



**FT880**



**FT922**