

# Gas Collection/Distribution

## Gas Distribution

### Gas Distribution Tubes, Pyrex

For dispersing gases in liquids. With sintered glass thimble on a tube stem. Thimble diameter 15mm. Overall length x tube diameter 300 x 8mm.

	Ref.	Thimble porosity
<b>GE700-10</b>	3830/00D	0
<b>GE700-11</b>	3830/02D	1
<b>GE700-12</b>	3830/04D	2
<b>GE700-13</b>	3830/06D	3

Absorption Tower – see *AB050*.


Gas Lighters, Flint and Piezo models – see *Burners section*.

Whatman In-Line filters for pipelines – see *Filtration section*.

### Pipeline Filter Tubes, Pyrex

With sintered glass flat disc set in-line with inlet and outlet tubes. Pressure differential across disc must not exceed 1 bar.

	Ref.	Disc porosity
<b>30mm disc, 10mm o.d. stem, 190mm long</b>		
<b>GE710-11</b>	3800/02M	1
<b>65mm disc, 12mm o.d. stem, 235mm long</b>		
<b>GE713-12</b>	3800/18D	2
<b>95mm disc, 14mm o.d. stem, 260mm long</b>		
<b>GE716-12</b>	3800/20M	2



GE700



GE710/713/716