

Microscopy

Slide Staining

Staining Block

Clear polished glass with flat glass cover. Cavity 6mm deep.

ML700-10 40 x 40mm

Staining Jar, Coplin pattern, glass

For holding 10 slides 76 x 26mm back to back. Complete with lid.

ML720-10 Glass

Staining Jar, Coplin pattern, plastic

Rugged polypropylene, with screw lid. For holding 10 slides 76 x 26mm back to back. Autoclavable but lid must be unscrewed and set gently on top.

ML720-55 Polypropylene

Staining Trough

Glass with loose flat lid. Grooved internally to hold 20 slides 76 x 26mm back to back.

ML730-10 Trough

Staining Trough

Glass with loose flat lid. Accepts one accessory tray. Overall 103 x 83 x 70mm W x D x H.

ML735-10 Trough

Accessory tray, glass. With inner grooves to accept 20 slides 76 x 26mm back to back.

ML737-08 Tray

Wire handle for ML737-08.

ML737-09 Wire handle

Slide rack ends

For suspending microscope slides over laboratory sinks using glass rod rails. Each end has two 7.5mm diameter holes with screw clamps to accommodate the glass rods, and two levelling screw feet. Supplied in pairs.

ML754-10 Staining rack ends

Slide Dispenser, Kartell

ABS plastic. Holds approximately 50 slides 76 x 26mm. Slides are ejected one at a time by turning either side knob.

ML798-15 Slide dispenser

Slide Drying Bench, Cole-Parmer SW200

With built in energy regulator and indicator lamps. Maximum temperature 100°C. Accepts up to fifty 76 x 26mm slides on removable drying bars. Overall 435 x 178 x 75mm W x D x H. For 220-240V a.c. single phase supplies, 150W. Formerly Electrothermal MH6616.

HG310-10 SW200


Histology equipment

Paraffin Section Mounting Bath – see HG300.

Paraffin Wax Dispenser – see HG320.




ML720-10, ML700, ML735, ML737, ML730



ML798 in use



ML720-55



ML754 in use



HG310 in use