

Pipettors

▲ SOCOREX SWISS Electronic Micropipettors

Acura Electro 926XS, 936, and 956 series

Lightweight, ergonomic, electronic pipettors with microprocessor control and motor driven to provide fatigue free, high performance pipetting.

- ◆ A choice of variable volume models up to 10,000µl and multichannel models up to 350µl
- ◆ Five programmable functions:
 - Pipetting
 - Step-by-step pipetting
 - Mixing
 - Reverse pipetting
 - Dilution
 - Tactile titration and gel loading
- ◆ Reversible display for left or right-handed use
- ◆ Three position pipetting speed selection
- ◆ Autoclavable and interchangeable pipette tip manifolds
- ◆ Adjustable tip ejector length accommodates various tip positions
- ◆ Easy user calibration mode

Acura 926XS and 936-series Single Channel Pipettors

As described. Supplied as a basic starter pack including pipettor, Q.C./calibration certificate, charging stand, a.c. adapter/charger and additional rechargeable battery, pipette tips and instructions.

	Range µl	Divisions µl	Accuracy % (at µl)	
Acura 926XS-series				
PR690-05	0.1 to 2.0	0.01	<±2.5 (0.1)	<±0.9 (2.0)
PR690-12	0.5 to 10	0.05	<±1.2 (0.5)	<±0.6 (10)
PR690-16	1.0 to 20	0.1	<±1.2 (1.0)	<±0.5 (20)
PR690-27	2.5 to 50	0.25	<±1.0 (2.5)	<±0.5 (50)
PR690-29	5 to 100	0.50	<±1.0 (5.0)	<±0.5 (100)
PR690-32	10 to 200	1.00	<±1.0 (5.0)	<±0.4 (200)
PR690-62	50 to 1000	5.00	<±0.8 (50)	<±0.4 (1000)
Acura 936-series				
PR690-83	100 to 2000	0.01	<±1.5 (100)	<±0.5 (2000)
PR690-87	250 to 5000	0.05	<±1.2 (250)	<±0.5 (5000)
PR690-92	500 to 10,000	0.05	<±1.0 (500)	<±0.5 (10,000)

Accessories

Pipettors only. With Q.C./calibration certificate, spare battery and instructions.

PR693-05	0.1 to 2.0µl
PR693-12	0.5 to 10µl
PR693-16	1.0 to 20µl
PR693-27	2.5 to 50µl
PR693-29	5 to 100µl
PR693-32	10 to 200µl
PR693-62	50 to 1000µl
PR693-83	100 to 2000µl
PR693-87	250 to 5000µl
PR693-92	500 to 10,000µl

Acura 956-series Multichannel Pipettors

As described. Supplied as a basic starter pack including pipettor, Q.C./calibration certificate, charging stand, a.c. adapter/charger and additional rechargeable battery, pipette tips and instructions.

	Range µl	Divisions µl	Accuracy % (at µl)	
Acura 956-series, 8-channel				
PR696-17	0.5 to 10	0.05	<±3.5 (0.5)	<±1.0 (10)
PR696-32	2.5 to 50	0.25	<±1.0 (2.5)	<±0.8 (50)
PR696-47	10 to 200	1.00	<±0.9 (10)	<±0.6 (200)
PR696-49	20 to 350	5.00	<±1.0 (20)	<±0.6 (350)
Acura 956-series, 12-channel				
PR696-52	0.5 to 10	0.05	<±3.5 (0.5)	<±1.0 (10)
PR696-67	2.5 to 50	0.25	<±1.0 (2.5)	<±0.8 (50)
PR696-82	10 to 200	1.00	<±0.9 (10)	<±0.6 (200)
PR696-90	20 to 350	5.00	<±1.0 (20)	<±0.6 (350)

Spares and accessories

PR694-15	Charging stand, 3-position (requires charger/a.c. adapter)
PR694-40	Charger/a.c. adapter 100-240V, with UK plug
PR694-45	Charger/a.c. adapter 100-240V, with European plug

Tips, polypropylene. Supplied in packs as indicated.

	Max. cap.	Colour	Pack qty.
PR580-09	10µl	Natural	1000
PR580-18	200µl	Yellow	1000
PR580-23	1000µl	Blue	1000
PR587-40	5000µl	Natural	250
PR587-55	10,000µl	Natural	100



PR690



PR690/696 with stand and battery



PR696-17 to PR696-49



PR696-52 to PR696-90