

Chromatography

Columns

With sintered disc porosity 0 and GP Rotaflo stopcock.

	<i>Dia.</i> <i>mm</i>	<i>Effective</i> <i>length</i> <i>mm</i>	<i>Socket</i> <i>size</i>	<i>GP Rotaflo</i> <i>stopcock</i> <i>bore, mm</i>
QCR10/20	10	200	14/23	3
QCR10/40	10	400	14/23	3
QCR20/30	20	300	19/26	3
QCR20/50	20	500	19/26	3
QCR30/50	25	500	24/29	3
QCR40/40	30	400	29/32	6
QCR60/50	40	500	40/38	6

Spare Rotaflo keys – see QGP series.



QCR10:-60: