

Extraction Liquid-Solid continued

Complete standard assemblies

	Nom. cap. ml	Extractor etc.	Components Flask	Condenser
Q100RASX	100	EX5/63	FR250/3S	CX7/06/SC
Q200RASX	200	EX5/83	FR500/3S	CX7/08/SC
Q600RDSX	600	EX5/75 DA5/75 JC75F	FR2L/5S	CX5/25/SC
Q2LRCSX	2000	EX5/105 DA5/100 JC100F	FR5L/5US	C13/55

Pressure equalising funnels

These funnels may be used with Soxhlets to facilitate recovery at the end of an extraction. Fitted between the extractor and the condenser the funnel stopcock is left open until extraction is completed. The stopcock is then closed so as to collect the solvent for future use. Full details see QDE series.

Paper thimbles

Recommended sizes for use with Quickfit extractors.

Extractor number	Recommended thimble size mm	Extractor number	Recommended thimble size mm
EX5/21	10 x 50	EX5/55/100	28 x 120
EX5/32	18 x 65	EX5/63	33 x 100
EX5/43	22 x 80	EX5/83	41 x 123
EX5/53	28 x 80	—	—
EX5/55	28 x 80	—	—

Note: The dimensions are internal and are given for single thickness thimbles only. If double thickness thimbles are used, the internal diameter should be 2mm less to allow for greater wall thickness.

Paper thimbles – see EX100/EX110.

Heating mantles – see EX200/EX300.

