

## Funnels

### Conical filter funnels

With sintered disc porosity 3. Semi-micro. Small scale for use with receivers with vents MD60 or MD61.

	Cone size	Disc dia. mm	Funnel dia. mm
<b>QSF4T31</b>	14/23	10	30

### Jointed receivers. With vent

	Capacity ml	Socket size	Height mm
<b>QMD60</b>	5	10/19	95

### Büchner filter funnels

With sintered disc porosity 3.

	Cone size	Disc dia. mm	Funnel dia. mm
<b>QSF3A32</b>	19/26	55	71
<b>QSF3A33</b>	24/29	55	71
<b>QSF3A34</b>	29/32	55	71
<b>QSF3B33</b>	24/29	95	106

### Powder funnels

Approximate funnel diameter: 105mm.

	Cone size
<b>QCF5/1</b>	14/23
<b>QCF5/2</b>	19/26
<b>QCF5/3</b>	24/29
<b>QCF5/4</b>	29/32
<b>QCF5/5</b>	34/35



QSF4T



QMD60



QMD61



QSF3



QCF5