# **Fractionation**

Socket size

19/26

19/26

## **Dufton columns**

QFC1/12

QFC1/22

These columns have an accurately fitted Dufton glass spiral which can be removed for cleaning purposes etc. This type of column has proved exceedingly successful in gas analysis.

Approx. effective

length

mm

150

300

Approx. overall

height

тт

270

425

Cone size

19/26

19/26

#### **Plain columns**

Designed for use with the various types of column packing which are available.

	Socket size	Cone size	Approx. effective length mm	Approx. overall height mm	
QFC7/23	24/29	24/29	250	365	
QFC7/43	24/29	24/29	500	615	

### Raschig ring packing for plain columns


Raschig rings are short lengths of glass tubing, usually having the same length as the diameter. Supplied in minimum quantities of 1 litre.

	Approx. diameter mm	Approx. length mm	
Raschig rings			
QFC8/09	9	9	

#### **Packing support**

These packing supports will hold the smallest column packing normally used.

	To fit socket size	Diameter mm	
QEX13/23	24/29	21	
QEX13/25	34/35	29	
QEX13/26	40/38	35	



QFC1:, QFC2:



QEX13:



0207 729 7998 e-mail: customer.services@thelabwarehouse.com fax orders: 020 7729 2657