

Stillheads

| | Flask | Cone size | Condenser |
|---------------------------------|-------|-----------|-----------|
| Recovery bends, sloping | | | |
| QSH1/11 | 14/23 | | 14/23 |
| QSH1/22 | 19/26 | | 19/26 |
| QSH1/23 | 24/29 | | 19/26 |
| QSH1/33 | 24/29 | | 24/29 |
| QSH1/44 | 29/32 | | 29/32 |
| Recovery bends, vertical | | | |
| QSH2/22 | 19/26 | | 19/26 |
| QSH2/33 | 24/29 | | 24/29 |

Plain

With thermometer sockets.

| | Flask | Cone size | Condenser | Socket size |
|---------|-------|-----------|-----------|-------------|
| QSH4/1 | 14/23 | | 14/23 | 14/23 |
| QSH4/22 | 19/26 | | 19/26 | 14/23 |
| QSH4/2 | 19/26 | | 19/26 | 19/26 |
| QSH4/23 | 24/29 | | 19/26 | 14/23 |
| QSH4/33 | 24/29 | | 24/29 | 14/23 |
| QSH4/44 | 29/32 | | 29/32 | 14/23 |

Claisen heads, sloping

| | | | | |
|---------|-------|--|-------|-------|
| QSH6/23 | 24/29 | | 19/26 | 14/23 |
| QSH6/33 | 24/29 | | 24/29 | 14/23 |
| QSH6/44 | 29/32 | | 29/32 | 14/23 |

Splash head adapter

| | Socket size | Cone size |
|----------|-------------|-----------|
| QSH17/23 | 19/26 | 24/29 |

Splash heads, vertical

| | Flask | Cone size | Condenser |
|---------|-------|-----------|-----------|
| QSH7/12 | 19/26 | | 19/26 |
| QSH7/13 | 24/29 | | 19/26 |
| QSH7/53 | 24/29 | | 24/29 |

Splash head adapters, rotary evaporator

For use with rotary evaporators to prevent "bumping" or foaming solvents from climbing to the condenser. Available with, or without, sintered glass (porosity 1) filter.

| | Flask | Cone size | Condenser |
|-----------------------|-------|-----------|-----------|
| Without filter | | | |
| QSH15/19 | 19/26 | | 29/32 |
| QSH15/24 | 24/29 | | 29/32 |
| QSH15/29 | 29/32 | | 29/32 |
| With filter | | | |
| QSH16/19 | 19/26 | | 29/32 |
| QSH16/29 | 29/32 | | 29/32 |



QSH1



QSH4



QSH17/23



QSH15



QSH2



QSH6



QSH7/11, QSH7/5



QSH16