

## Stillheads

	Flask	Cone size	Condenser
<b>Recovery bends, sloping</b>			
QSH1/11	14/23		14/23
QSH1/22	19/26		19/26
QSH1/23	24/29		19/26
QSH1/33	24/29		24/29
QSH1/44	29/32		29/32
<b>Recovery bends, vertical</b>			
QSH2/22	19/26		19/26
QSH2/33	24/29		24/29

### Plain

With thermometer sockets.

	Flask	Cone size	Condenser	Socket size
QSH4/1	14/23		14/23	14/23
QSH4/22	19/26		19/26	14/23
QSH4/2	19/26		19/26	19/26
QSH4/23	24/29		19/26	14/23
QSH4/33	24/29		24/29	14/23
QSH4/44	29/32		29/32	14/23

### Claisen head, sloping

QSH6/44	29/32		29/32	14/23
---------	-------	--	-------	-------

### Splash head adapter

	Socket size	Cone size
QSH17/23	19/26	24/29

### Splash heads, vertical

	Flask	Cone size	Condenser
QSH7/12	19/26		19/26
QSH7/13	24/29		19/26
QSH7/53	24/29		24/29

### Splash head adapters, rotary evaporator

For use with rotary evaporators to prevent "bumping" or foaming solvents from climbing to the condenser. Available with, or without, sintered glass (porosity 1) filter.

	Flask	Cone size	Condenser
<b>Without filter</b>			
QSH15/19	19/26		29/32
QSH15/24	24/29		29/32
QSH15/29	29/32		29/32
<b>With filter</b>			
QSH16/19	19/26		29/32
QSH16/29	29/32		29/32



QSH1



QSH4



QSH17/23



QSH15



QSH2



QSH6/44



QSH7/12, QSH7/13



QSH16