

Refrigeration

Pharmacy Refrigerators



For the storage of medicines and vaccines under controlled conditions in dispensaries, surgeries and in general use. Maintain internal air temperature between +2°C and +8°C.

- ◆ Intellicold™ touchscreen controller ensures accuracy and consistency within the specified temperature range whilst only using the power necessary to do so, giving substantially lower running costs
- ◆ Built-in temperature display with max./min. recordings and continuous readout, high/low temperature alarm with reset, power failure and "door open" alarms
- ◆ Fan assisted air circulation provides enhanced uniformity
- ◆ Conform with UK Dept. of Health "Green Book", Royal Pharmaceutical Society (RPS), General Pharmaceutical Council (GPHC) and National Travel Health Network and Centre (NaTHNaC) guidelines
- ◆ Integral door lock and internal lighting (glass door models have independent external light switch)
- ◆ Remote alarm contacts
- ◆ 13 month temperature data storage with expansion via SD card slot (accessory coded SD card available)
- ◆ R134a, CFC-, HCFC- and ammonia-free refrigerant and environmentally friendly insulation used
- ◆ Moulded plastic interior with enamelled steel exterior

Pharmacy Refrigerators

As described. With automatic defrost, lockable door, Intellicold™ touchscreen controller, fan air circulation, wire shelves and solid or glass-paned doors as indicated. For 230V 50Hz single phase supplies.

	Gross cap. litres	H x W x D mm	Weight kg	Shelves supplied
Solid door models				
RG270-17	36	540x 450 x 510	30	2
RG270-29	66	735 x 450 x 510	34	3
RG270-37	150	835 x 600 x 600	51	3
RG270-47	340	1500 x 600 x 700	79	6
RG270-52	440	1865 x 600 x 700	93	6
RG270-71	543	1960 x 750 x 755	100	5
Glass door models				
RG270-72	150	835 x 600 x 600	54	3
RG270-89	340	1500 x 600 x 700	91	6
RG270-92	440	1865 x 600 x 700	105	6
RG270-96	543	1960 x 750 x 755	100	5

RG270-99 Accessory coded SD card for RG270-series only, allows up to 13 months data storage

Note:

Due to the large size and weight of these units the user must ensure that suitable space area, floor loading and access for installation is available prior to ordering. They are installable by the user and delivery therefore will only be made to the customers unloading area.



RG270-37 in use



RG270-72 in use



☎ 0207 729 7998

e-mail: customer.services@thelabwarehouse.com fax orders: 020 7729 2657

