Waste, Watch Glasses

Beta and Gamma Radiation Shielding, Boxes

A range of benchtop storage boxes with hinged lids to accept accessory inserts as indicated. Dimensions stated are H x W x D.

For beta radiation

RA612-15	Mini-box, 75 x 105 x 105mm					
RA612-25	Midi-box,					
	80 x 185 x 105mm					
RA612-35	Maxi-box,					
	160 x 300 x 185mm					
For gamma radiation						
RA615-25	Midi-box,					
	784 x 189 x 109mm					

RA615-25 Midi-box, 784 x 189 x 109mm RA615-35 Maxi-box, 164 x 204 x 180mm

164 x 304 x 189mm

SAFETY NOTES

1) Beta radiation shields are unsuitable for use against secondary x-rays (Bremsstrahlung) or gamma radiation.

2) Gamma radiation shields will block emissions effectively from ¹²⁵I and any low energy gamma emitters. They are not suitable for more energetic isotopes of iodine. Heavierwalled shields are available for use with these isotopes - details on request. Gamma shields must not be used with beta isotopes since Bremsstrahlung will be produced.

3) Radioactive waste must never be stored on the benchtop for long periods and must be disposed of properly. Proper handling and storage of isotopes minimises unnecessary exposure to radiation.

Watch Glasses

Watch Glasses

Clear glass with ground edges. Diameters are nominal. Supplied singly.

	Dia.	
	mm	
WG100-10	30	
WG100-15	40	
WG100-20	50	
WG100-25	60	
WG100-30	70	
WG100-35	80	
WG100-40	90	
WG100-45	100	
WG100-50	130	
WG100-55	150	
WG100-60	200	

Watch glasses

Polypropylene. Transparent with stable ring base. Autoclable 121°C.

	Dia. mm
WG145-23	60
WG145-38	80
WG145-48	100
WG145-52	125

Gas jar covers - see GE302.

Weighing boats - see BD750/BD757.


Watch Glasses, Pyrex

Borosilicate glass. Diameters are nominal. Supplied

singly.					
	Ref. 1060-	Dia. mm			
WG130-13	/04D	40			
WG130-18	/06D	50			
WG130-23	/09D	60			
WG130-38	/10D	80			
WG130-48	/14D	100			
WG130-52	/18D	125			
WG130-58	/20D	150			
WG130-63	/21D	200			
WG130-70	/24D	250			



RA612-15 in use with accessories



RA612-25 with accessories, RA615-25 is similar











© 0808 200 5227 e-mail: sales@thelabwarehouse.com fax orders: 020 7729 2657

