

# Water Testing

## Water Sample Bottles



Clear PET, with choice of 350ml, 500ml or 1000ml capacities and stable, square or space-saving, rectangular profile. With leakproof, 'double cap' with tamper-evident seal and colour-coded blue (dosed with 120mg/l sodium thiosulphate) or white (undosed). Both versions are labelled with batch details and expiry date. Sterilised by gamma irradiation. Supplied in packs as indicated.

	Ref.	Capacity ml	Dosed?	Cap	Pack qty.
<b>Rectangular</b>					
<b>WR284-02</b>	WSC350	350	✓	Blue	144
<b>WR284-04</b>	WSC350/NT	350	—	White	144
<b>WR284-06</b>	WSC500T	500	✓	Blue	144
<b>WR284-08</b>	WSC500T/NT	500	—	White	144
<b>Square</b>					
<b>WR284-32</b>	WSC500	500	✓	Blue	108
<b>WR284-34</b>	WSC500/NT	500	—	White	108
<b>WR284-36</b>	WSC1000	1000	✓	Blue	70
<b>WR284-38</b>	WSC1000/NT	1000	—	White	70

## Water Sampling Bottles

Sterile, rectangular, HDPE bottles in a range of sizes with offset neck for easy dispensing, tamper-evident caps and attached label. Supplied dosed with Sodium Thiosulphate or undosed and in packs as indicated.

	Capacity ml	Dosed?	Cap	Pack qty.
<b>WR286-10</b>	250	✓	Blue	140
<b>WR286-14</b>	250	—	Green	140
<b>WR286-20</b>	500	✓	Blue	84
<b>WR286-24</b>	500	—	Green	84
<b>WR286-30</b>	1000	✓	Blue	56
<b>WR286-34</b>	1000	—	Green	56



WR284



WR286

## Thio-Bags®



Whirl-Pak® sample bags suitable for water sampling. Available in either non self-standing or self-standing versions. Ready dosed with sodium thiosulphate to neutralise chlorine in the sample (100ml capacity bags contain 10mg and 300ml capacity bags 30mg of active sodium thiosulphate). Supplied in packs of 100.

### Non Self-standing

**WR290-02** Thio-Bags®, 100ml

**WR290-11** Thio-Bags®, 300ml

### Self-standing

**WR290-07** Thio-Bags®, 100ml

**WR290-12** Thio-Bags®, 300ml

### Racks for Thio-Bags®

Dimensions shown are number of apertures x dimensions of apertures.

**BA110-10** 21 x (40 x 50mm), for WR290-02/-07

**BA110-20** 18 x (50 x 50mm), for WR290-02/-07

**BA110-30** 20 x (64 x 64mm), for WR290-11/-12



WR290-02, -11 is similar



BA110-10 in use



BA110-30 in use



WR290-07, -12 is similar