

Autoclave Accessories

Autoclave Accessories

Autoclave Deodoriser Capsules

Eliminate offensive odours generated from autoclaving organic waste. Available in a choice of fragrances as indicated. One or two capsules is sufficient to deodorise one load. Supplied in packs of 100 capsules.

- AU604-35** Rose fragrance
- AU604-40** Citrus fragrance

Autoclave Tape

Self adhesive paper with indicating stripes which change colour after being autoclaved. Conforms to European Standard EN867. Dimension stated is tape width. Supplied in rolls of 55 metres.

- AU630-10** 18mm
- AU630-12** 24mm
- AU630-95** Tape dispenser for 18/24mm wide autoclave tape

Sterilising Pouches, self-seal/transparent

Equally suitable for autoclave or ethylene oxide gas sterilisation. The pouch has a two-ply polyester/ polypropylene laminated transparent front bonded to a medical grade white sterilisation paper. The paper used conforms to BS6256 and the complete pouch meets the requirements of DHSS specification TSS/S/ 330.010. Dual arrows are printed on each pouch; pink changes to brown under normal autoclaving conditions and blue changes to yellow when gas sterilised.

A self-seal adhesive band, which is on the shortest side, is used to seal the pouch. Heat sealing is not required.

	Pouch size mm	Pack qty
AU660-17	90 x 203	200
AU660-22	90 x 230	200
AU660-32	130 x 254	200
AU660-37	130 x 360	200
AU660-42	190 x 330	200
AU660-47	300 x 390	200

Biohazard bags suitable for autoclaving – see *Bags* section.



AU604



AU630

Sterilising Bags, paper

For autoclave sterilising. To BS6257 and using medical grade white sterilisation paper to BS6256. Available with alternative closures:-

Heat seal (AU690 series) – supplied with an adhesive band on the inside of the bag for sealing using a suitable heat seal machine. A pink indicator on each bag changes to brown under normal autoclaving conditions.

Plain closure (AU695 series) – are folded and secured using the accessory indicating autoclave tape.

All bags are gusseted and strongly constructed with a double turn at the bottom for extra safety. The opening is on the shortest side.

W x L mm	Gusset mm	Pack qty	Heat seal closure	Plain top closure
90 x 170	50	1000	AU690-15	AU695-16
110 x 190	30	1000	—	AU695-21
140 x 250	75	1000	AU690-52	AU695-52
140 x 330	50	500	AU690-55	AU695-55
180 x 380	95	500	AU690-64	AU695-64
180 x 530	75	250	—	AU695-67
190 x 330	65	500	AU690-74	—
250 x 380	100	250	—	AU695-79
300 x 610	125	250	—	AU695-88

Autoclave tape for securing the plain top closure bags – see *AU630*.

Heat resistant gloves – see *SA704* in the *Safety* section.

Sealing film – see *SE162*.



SA704



AU660

AU690, AU695