Sample Bags

WHIRL-PAK°

Polyethylene with tab and wire frame at top which allows bag to be opened wide. The bag is sealed by pulling the tabs and whirling the bag around the wire. Can be opened and closed as required. Ethylene oxide gas sterilised. Suitable for solids and liquids. Available with or without writing panel. Supplied in

D x W,	Capacity,	Plain	Writing	Writing panel,	Pack
mm	ml		Panel	Self-standing	qty.
125 x 65	29		BA105-11		500
125 x 75	58	ВА104-13	BA105-13		500
185 x 75	118	BA104-15	BA105-15	BA107-15	500
185 x 75	118	—	BA106-15*	—	500
180 x 95	207	BA104-16	BA105-16	—	500
230 x 115	532	BA104-17	BA105-17	BA107-17	500
230 x 115	532	—	BA106-17*	—	500
230 x 150	710	BA104-19	BA105-19	BA107-19	500
230 x 150	710	BA104-21**	_	BA107-25 [†]	500
305 x 125	798	BA104-23	_		500
380 x 125	1065	BA104-25	_		500
380 x 150	1242	BA104-27	BA105-29	BA107-27 [†]	500
380 x 190	2041	BA104-29		BA107-29 [†]	500
380 x 254	2721		BA105-31	_	250
508 x 254	3637		BA105-33	_	250
508 x 380	5441		BA105-36	_	100

^{*}With flat wire band.

Additional Whirl-Pak® bags:

Hygiene sampling bags - see HT845-/HT850-series.

Blender bags - see MS865-series.

Water sample test bags - see WR290-series.

Locking Pipe Closures, Nasco

Autoclavable plastic, locking tubes that allow temporary air- and water-tight closure of Whirl-Pak bags without whirling the bag and folding the tabs. The bag is wrapped around the inner core and clamped using the outer sleeve. Can also be used with other sample bags up to 190mm wide. Supplied in a pack of 50.

BA110-08 Locking pipe closures

Whirl-Pak® Wire racks, Nasco

For convenient storage or transportation of samples within Whirl-Pak® bags. Dimensions shown are number of apertures x dimensions and to hold bags with capacities as indicated.

BA110-10 21 x (40 x 50mm), up to 207ml bags **BA110-20** 18 x (50 x 50mm), up to 207ml bags **BA110-30** 20 x (64 x 64mm), 532-1065ml bags **BA110-40** 6 x (114 x 114mm), 710-2041ml bags **BA110-60** Carrying handle for BA110-racks

Whirl-Pak® Scoop bags, Nasco

Comprises 532ml Whirl-Pak® bag, 230 x 115mm, containing a sterile, polystyrene, 60ml scoop. Supplied in a pack of 50.

BA115-10 Scoop bags

Biohazard Specimen bags, Nasco

Whirl-Pak® Trans-Pak® clinical specimen transport bags, non-sterile. Bag dimensions 230 x 150mm. With external pouch 146 x 137mm and biohazard label. Accessory security tape is made of destructible material which cannot be removed in one piece and immediately tears if removal is attempted, indicating any tampering. Supplied in packs as indicated.

BA118-15 Biohazard specimen bags, pk.500 **BA118-17** Security tape strips, pk.100

Soil Sample bags, Nasco

Used in sample submission for analytical assessment. Bleached paper with polyethylene liner and heat sealed base so that no sample sifting occurs. Capacity 769ml. With attached, folding tin strap tie. Supplied in a pack of 500.

BA120-25 Soil sample bags







BA105 in use, BA104 are similar



BA107-27 in use



BA115-08





BA118-15



BA118-17



BA120-25 in use



C +44 (0)20 7729 2405

BA110-40 in use with BA110-60

^{**}With separate external pocket for holding sample paperwork.

[†]Supplied in a pack of 250.