

# Bottles, Glass, Capped

## Glass Bottles with Caps continued

### Amber winchester pattern

With black polyethylene cap and cone liner.

	Cap. ml	Neck/cap	Pack qty.
<b>BR372-50</b>	500	R6/31	24
<b>BR372-80</b>	1000	R6/31	12
<b>BR372-90</b>	2000	R6/31	12
<b>BR372-95</b>	2500	R6/31	12

### Wide mouth powder

With black plastic cap and PVDC liner.

	Cap. ml	Neck/cap	Pack qty.
<b>Clear glass</b>			
<b>BR400-16</b>	30ml	R3/33	80
<b>BR400-24</b>	60ml	R3/38	72
<b>BR400-33</b>	4oz	R3/38	80
<b>BR400-36</b>	175ml	R3/48	53
<b>BR400-41</b>	250ml	R3/48	42
<b>BR400-49</b>	500ml	R3/58	33

### Amber glass

<b>BR402-17</b>	30ml	R3/33	127
<b>BR402-23</b>	60ml	R3/38	72
<b>BR402-34</b>	4oz	R3/38	68
<b>BR402-41</b>	250ml	R3/48	42
<b>BR402-50</b>	500ml	R3/58	28
<b>BR402-80</b>	1000ml	R3/63	12

### Wide mouth powder

With white polypropylene wadless cap.

	Cap. ml	Neck/cap	Pack qty.
<b>Clear glass</b>			
<b>BR404-33</b>	4oz	R3/38	80
<b>BR404-41</b>	250	R3/48	42
<b>Amber glass</b>			
<b>BR406-34</b>	4oz	R3/38	68
<b>BR406-41</b>	250ml	R3/48	42

### Ointment jar

Extra wide mouth with black plastic cap.

	Cap. ml	Neck/cap	Pack qty.
<b>Clear glass</b>			
<b>BR410-24</b>	60	R3/51	72
<b>BR410-30</b>	120	R3/58	54
<b>BR410-42</b>	250	R3/83	48
<b>BR410-50</b>	500	R3/100	24

### Important Warning Sterilising Bottles

Some bottles where indicated can be sterilised by autoclaving. However BEFORE AUTOCLAVING (or sterilising by gas or dry heat) the cap or closure must either be removed or just set on top of the container **WITHOUT ENGAGING THE THREAD.**

Alternative bottle/cap combinations are available – details on request.



BR372



BR402



BR406



BR400



BR404



BR410