

Centrifuge Tubes

Nalgene Tubes and Bottles

All dimensions and capacities are nominal.

Conical, polypropylene, graduated, rimless

Autoclavable. For non-refrigerated low speeds up to 6000xg. Supplied in packs of 10.

| | Ref. 3103- | Cap. ml | L x dia. mm |
|----------|---------------|------------|----------------|
| CF604-22 | 0015 | 15 | 119 x 17 |
| CF604-32 | 0050 | 50 | 134 x 28.5 |

Conical, polycarbonate, graduated, rimless

Transparent, autoclavable. For low speeds up to 6000xg. Supplied in packs of 10.

| | Ref. 3105- | Cap. ml | L x dia. mm |
|----------|---------------|------------|----------------|
| CF608-18 | 0015 | 15 | 120.2 x 17 |
| CF608-28 | 0050 | 50 | 135 x 29 |

Straight, polypropylene, rimless

Autoclavable. For non-refrigerated speeds up to 50000xg. Supplied in packs of 10.

| | Ref. 3110- | Cap. ml | L x dia. mm |
|----------|---------------|------------|----------------|
| CF640-16 | 0120 | 12 | 103 x 16 |
| CF640-18 | 0150 | 15 | 114 x 16 |
| CF640-20 | 0160 | 16 | 100 x 18 |
| CF640-34 | 0500 | 50 | 103 x 28.7 |

Straight, polycarbonate, rimless

Transparent, autoclavable. For use at speeds up to 50000xg. Supplied in packs of 10.

| | Ref. 3117- | Cap. ml | L x dia. mm |
|----------|---------------|------------|----------------|
| CF655-16 | 0120 | 12 | 104 x 16 |
| CF655-18 | 0150 | 15 | 114.5 x 16 |
| CF655-20 | 0160 | 16 | 101 x 18 |
| CF655-34 | 0500 | 50 | 104 x 29.5 |

Bottles, polycarbonate

Transparent, autoclavable. With polypropylene narrow mouth screw cap. Supplied in packs of 4.

Maximum speed ratings:

| | | |
|--------|---|---------|
| 250ml | - | 27500xg |
| 500ml | - | 13700xg |
| 1000ml | - | 7100xg |

| | Ref. 3122- | Cap. ml | L x dia. mm |
|----------|---------------|------------|----------------|
| CF830-20 | 0250 | 250 | 128 x 62 |
| CF830-30 | 0500 | 500 | 170 x 69.5 |
| CF830-40 | 1000 | 1000 | 188.4 x 97.6 |

Important Warning Sterilising Centrifuge Ware

Some tubes and bottles where indicated can be sterilised by autoclaving. However BEFORE AUTOCLAVING (or sterilising by gas or dry heat) the cap or closure must either be removed or just set on top of the container WITHOUT ENGAGING THE THREAD.

Nalgene Oak Ridge Tubes

All dimensions and capacities are nominal.

All tubes are autoclavable and can be used in refrigerated and non-refrigerated centrifuges up to 50000xg.

Supplied with non-contaminating screw cap without liner which is leakproof during ordinary use outside the centrifuge.

Length x diameter:

| | |
|------|------------|
| 10ml | 81 x 16mm |
| 30ml | 95 x 26mm |
| 50ml | 107 x 29mm |

Polycarbonate, transparent, with polypropylene screw cap


| | Ref. 3118- | Cap. ml | Pack qty. |
|----------|---------------|------------|--------------|
| CF900-10 | 0010 | 10 | 10 |
| CF900-30 | 0030 | 30 | 10 |
| CF900-45 | 0050 | 50 | 10 |

Polyallomer, transparent, with polypropylene screw cap

| | Ref. 3119- | Cap. ml | Pack qty. |
|----------|---------------|------------|--------------|
| CF915-10 | 0010 | 10 | 10 |
| CF915-30 | 0030 | 30 | 10 |
| CF915-45 | 0050 | 50 | 10 |

PTFE/FEP, with ETFE screw cap

| | Ref. 3114- | Cap. ml | Pack qty. |
|----------|---------------|------------|--------------|
| CF935-10 | 0010 | 10 | 2 |
| CF935-18 | 0030 | 30 | 2 |
| CF935-25 | 0050 | 50 | 2 |



CF604



CF608



CF655



CF900



CF640



CF830



CF915, CF935 is similar