Cleaning

Disinfectants

Microbicidal Surface Disinfectant, Terralin® Protect

High-performance, broad spectrum disinfectant and cleaner, effective against bacteria (including Tb bacilli), fungi, MRSA, Hepatitis B viruses, rota viruses, HIV and adeno viruses. For use on large surfaces, fittings and equipment that can be wet-wiped.

- Non-toxic, low-odour formulation
- Economical low concentration required for effective use (0.25% to 2% in water depending on application/exposure required)
- Good cleaning and soil-loosening properties

As described. Supplied in pack as indicated.

CL850-10 Terralin® Protect, 2 litres

Microbicidal Equipment Disinfectant, Mikrozid® AF

Alcohol-based, rapid cleaner and disinfectant for medical devices and equipment in clinical areas. Effective against Bacteria (including Tb bacilli), fungi, hepatitis B viruses, HIV, HCV, adeno viruses, MRSA and Polio.

- Easy-to-use
- Aldehyde-free
- Leaves no residue

A single application and thorough wiping with a disposable tissue will ensure effective cleaning. Secondary application and leaving the surface to dry will provide disinfection.

Note: Mikrozid should not be applied to alcohol-sensitive surfaces such as Perspex (acrylic). Supplied in pack as indicated.

CL860-15 Mikrozid AF, 1 litre

CL862-04 Spray nozzle, for use with CL860-15





CL860-15 with CL862-04





