

Flasks, Glass

Flask Supports

Cork Rings, compressed

For supporting round bodied flasks. Depth as indicated. Supplied singly.

	<i>o.d. x i.d. mm</i>	<i>For flasks ml</i>
35mm deep		
CW308-08	80 x 30	25 to 150
CW308-12	110 x 60	250, 500
CW308-16	140 x 90	1000
CW308-20	170 x 120	2000, 3000
CW308-24	210 x 150	4000, 5000
45mm deep		
CW309-30	240 x 200	4000 to 10000

Flask Support Stand

Polypropylene. For round bottom flasks. Top diameter 145mm.

FK180-15 Flask support stand

Flask Support Rings

Blue silicone rubber. Operating range -60 to +200°C. Double sided. Sit firmly on the bench and with bowl either side to hold different size flasks. Hold flask firmly but gently, irrespective of the angle.

Size	Bowl diameters	
	Wide	Narrow
1	76mm	44mm
2	121mm	80mm
3	165mm	118mm

Supplied as a set of 3 (one of each size)

FK180-60 Support rings set

Spare support rings for FK180-60

Supplied singly.

FK180-62 Support ring 76/44mm

FK180-64 Support ring 121/80mm

FK180-66 Support ring 165/118mm

Flask Stabilising Rings

Vinyl-coated lead. Wrap around flasks up to 4000ml to provide stability whilst on a bench or immersed in a water bath. Choice of closed, or open ring format which allows slight adjustment of the ring opening and position. An accessory holder is available to hold up to 8 lead rings. For flask sizes as indicated. Supplied singly.

	<i>Flask sizes, ml</i>	<i>i.d. mm</i>	<i>o.d. mm</i>	<i>Weight, g</i>
Closed Ring				
FK190-10	125 to 500	48	77	277
FK190-12	250 to 1000	51	90	508
FK190-14	500 to 2000	57	98	635
FK190-16	1000 to 4000	70	115	970
Open Ring				
FK192-10	125 to 500	42	75	227
FK192-12	250 to 1000	51	90	472
FK192-14	500 to 2000	54	98	599
FK192-16	1000 to 4000	60	113	907



CW308/CW309



FK180-15



FK190/FK192 in use



FK180-60 in use