

# Flasks, Glass

## Flask Supports

### Cork Rings, compressed

For supporting round bodied flasks. Depth as indicated. Supplied singly.

	<i>o.d. x i.d. mm</i>	<i>For flasks ml</i>
<b>35mm deep</b>		
<b>CW308-08</b>	80 x 30	25 to 150
<b>CW308-12</b>	110 x 60	250, 500
<b>CW308-16</b>	140 x 90	1000
<b>CW308-20</b>	170 x 120	2000, 3000
<b>CW308-24</b>	210 x 150	4000, 5000
<b>45mm deep</b>		
<b>CW309-30</b>	240 x 200	4000 to 10000

### Flask Support Stand

Polypropylene. For round bottom flasks. Top diameter 145mm.

**FK180-15** Flask support stand

### Flask Support Rings

Blue silicone rubber. Operating range -60 to +200°C. Double sided. Sit firmly on the bench and with bowl either side to hold different size flasks. Hold flask firmly but gently, irrespective of the angle.

Size	Bowl diameters	
	Wide	Narrow
1	76mm	44mm
2	121mm	80mm
3	165mm	118mm

Supplied as a set of 3 (one of each size)

**FK180-60** Support rings set

### Spare support rings for FK180-60

Supplied singly.

**FK180-62** Support ring 76/44mm

**FK180-64** Support ring 121/80mm

**FK180-66** Support ring 165/118mm

## Flask Stabilising Rings

Vinyl-coated lead. Wrap around flasks up to 4000ml to provide stability whilst on a bench or immersed in a water bath. Choice of closed, or open ring format which allows slight adjustment of the ring opening and position. An accessory holder is available to hold up to 8 lead rings. For flask sizes as indicated. Supplied singly.

	<i>Flask sizes, ml</i>	<i>i.d. mm</i>	<i>o.d. mm</i>	<i>Weight, g</i>
<b>Closed Ring</b>				
<b>FK190-10</b>	125 to 500	48	77	277
<b>FK190-12</b>	250 to 1000	51	90	508
<b>FK190-14</b>	500 to 2000	57	98	635
<b>FK190-16</b>	1000 to 4000	70	115	970
<b>Open Ring</b>				
<b>FK192-10</b>	125 to 500	42	75	227
<b>FK192-12</b>	250 to 1000	51	90	472
<b>FK192-14</b>	500 to 2000	54	98	599
<b>FK192-16</b>	1000 to 4000	60	113	907



**CW308/CW309**



**FK180-15**



**FK190/FK192 in use**



**FK180-60 in use**