

# Pipettes, One-Mark

## Glass, One- and Two-Mark to ISO648 & DIN12690

- ◆ Colour coded to BS3996
- ◆ Calibrated for delivery
- ◆ Class AS pipettes only require a 5 second waiting time
- ◆ 1ml and 2ml sizes are straight pattern without bulb
- ◆ 15ml sizes are supplementary to the DIN12690 and ISO648 ranges

Colour code					
1ml	Blue	10ml	Red	25ml	Blue
2ml	Orange	15ml	Green	50ml	Red
5ml	White	20ml	Yellow	100ml	Yellow

Tolerances		
Cap. ml	Class B $\pm ml$	Class AS $\pm ml$
1	0.015	0.008
2	0.020	0.010
5	0.030	0.015
10	0.040	0.020
15	0.050	0.025
20	0.060	0.030
25	0.060	0.030
50	0.100	0.050
100	0.150	0.080

## MBL®

### One-Mark Pipettes, MBL

Soda-lime glass. Amber stain graduations and inscriptions.

	Ref.	Cap. ml
<b>Class B</b>		
<b>PM102-10</b>	PRA700B	1
<b>PM102-13</b>	PRA702B	2
<b>PM102-22</b>	PRA704B	5
<b>PM102-25</b>	PRA706B	10
<b>PM102-30</b>	PRA708B	15
<b>PM102-35</b>	PRA710B	20
<b>PM102-40</b>	PRA712B	25
<b>PM102-45</b>	PRA714B	50
<b>PM102-50</b>	PRA716B	100
<b>Class AS</b>		
<b>PM107-10</b>	PRA720AS	1
<b>PM107-13</b>	PRA722AS	2
<b>PM107-22</b>	PRA724AS	5
<b>PM107-25</b>	PRA726AS	10
<b>PM107-30</b>	PRA728AS	15
<b>PM107-35</b>	PRA730AS	20
<b>PM107-40</b>	PRA732AS	25
<b>PM107-45</b>	PRA734AS	50
<b>PM107-50</b>	PRA736AS	100
<b>Class AS Works Certified</b>		
<b>PM108-10</b>	PRA740ASWC	1
<b>PM108-13</b>	PRA742ASWC	2
<b>PM108-22</b>	PRA744ASWC	5
<b>PM108-25</b>	PRA746ASWC	10
<b>PM108-30</b>	PRA748ASWC	15
<b>PM108-35</b>	PRA750ASWC	20
<b>PM108-40</b>	PRA752ASWC	25
<b>PM108-45</b>	PRA754ASWC	50
<b>PM108-50</b>	PRA756ASWC	100


### Two-Mark Pipettes, MBL

Soda-lime glass. Amber stain graduations and inscriptions.

	Ref.	Cap. ml
<b>Class AS</b>		
<b>PM115-10</b>	PRA620AS	1
<b>PM115-13</b>	PRA622AS	2
<b>PM115-22</b>	PRA624AS	5
<b>PM115-25</b>	PRA626AS	10
<b>PM115-35</b>	PRA630AS	20
<b>PM115-40</b>	PRA632AS	25
<b>PM115-45</b>	PRA634AS	50
<b>PM115-50</b>	PRA636AS	100



PM102



PM107, PM108/PM115 are similar