

Pipettors

Microlitre Transferpette®



- ◆ Ergonomically designed
- ◆ Insulating jacket for increased precision
- ◆ Choice of models
 - digital variable volumes 0.1µl to 5000µl
 - fixed volumes 5µl to 2000µl
- ◆ Conformity certified
- ◆ Integral tip ejector

The lightweight construction and ergonomic design make it possible to hold your hand quite naturally when pipetting. This helps to overcome the fatigue and strain of repetitive pipetting.

The insulating jacket largely protects the measuring unit against external temperature effects, such as the warmth of the hand.

The controls are arranged so that they cannot be confused, even when working at speed. The tip ejection key is different from the pipetting key in both position and appearance. The compact casing has neither sharp edges nor a rough surface. Working with the Transferpette® is comfortable, even when wearing gloves.

Transferpette®, variable volume models

Digital reading.

	Range µl	Subdiv. µl	A%±	V%±	Tips reqd.
PR134-06	0.1 to 1.0	0.005	2.0	1.2	Nano-cap
PR134-08	0.5 to 10	0.05	1.0	0.8	Crystal
PR134-12	2 to 20	0.10	0.8	0.4	Crystal
PR134-17	2 to 20	0.10	0.8	0.4	Yellow
PR134-25	5 to 50	0.1	0.8	0.4	Yellow
PR134-30	10 to 100	0.1	0.6	0.2	Yellow
PR134-32	20 to 200	1.0	0.6	0.2	Yellow
PR134-35	25 to 250	1.0	0.6	0.2	Blue
PR134-40	100 to 1000	1.0	0.6	0.2	Blue
PR134-45	500 to 5000	10.0	0.6	0.2	Natural

Transferpette®, fixed volume models

	Volume µl	A%±	V%±	Tips reqd.
PR138-22	5	1.0	0.8	Yellow
PR138-25	10	1.0	0.8	Yellow
PR138-28	20	0.8	0.4	Yellow
PR138-32	25	0.8	0.4	Yellow
PR138-35	50	0.8	0.4	Yellow
PR138-38	100	0.6	0.2	Yellow
PR138-40	200*	0.6	0.2	Yellow
PR138-42	200*	0.6	0.2	Blue
PR138-45	250	0.6	0.2	Blue
PR138-48	500	0.6	0.2	Blue
PR138-52	1000	0.6	0.2	Blue
PR138-56	2000	0.6	0.2	Natural

*Two variants of the popular 200µl volume models are available depending on the tip style preferred.

A% = Accuracy } Figures quoted are based on the maximum
V% = Reproducibility } volume and distilled water at 20°C.

Accessory Racks

- PR142-08 Bench mounting for up to six Transferpettes.
Divides in half for use against an upright surface
- PR142-10 Wall mounting for up to three Transferpettes

Disposable Pipette Tips, Brand

High quality polypropylene tips which are Conformity Certified when used with Transferpette® pipettors. Cadmium free. Also suitable for use with Eppendorf and Gilson pipettors.

TipBox Set

Supplied as a steam sterilisable plastic box with hinged transparent cover containing a filled tip rack, plus additional tip racks as indicated.

TipRacks

Refill packs for TipBox Sets, comprising racks for tips to be dropped into empty TipBoxes as required.

Loose Packs

Supplied loose in a bag or carton in the quantities indicated.

Nano-Cap™, for volumes 0.1 to 20µl

- PR140-05 TipBox set, 5 racks x 96 tips
- PR140-08 TipRacks for PR140-05, 10 racks x 96 tips
- PR140-10 Loose pack of 2000 tips

Crystal, for volumes 0.5 to 20µl

- PR140-12 TipBox set, 5 racks x 96 tips
- PR140-14 TipRacks for PR140-12, 10 racks x 96 tips
- PR140-16 Loose pack of 2000 tips

Yellow, for volumes 2 to 200µl

- PR140-24 TipBox set, 5 racks x 96 tips
- PR140-26 TipRacks for PR140-24, 10 racks x 96 tips
- PR140-28 Loose pack of 1000 tips

Natural*, for use with PR134-08 and -12 only

- PR140-45 TipBox set, 5 racks x 96 tips
- PR140-48 TipRacks for PR140-45, 10 racks x 96 tips
- PR140-52 Loose pack of 1000 tips

*Except for PR140-52 which are blue coloured tips.

Natural, for use with PR134-45 and PR138-56 only

- PR140-80 Loose pack of 200 tips




PR134, PR138 is similar with fixed volume



PR142-08 and PR142-10 with instruments

Conformity Certified

Indicated by the symbol  shows that the product has been manufactured by Brand under quality control procedures according to Eichordnung the German Federal Weights and Measures Regulations.

Provided the instrument is used with the recommended accessories and in accordance with the supplied instructions, Brand certify conformity with annex 12 of the Eichordnung.