

# Pipettors

## Multichannel Transferpette® S - 8/ -12



- ◆ Conformity certified
- ◆ For pipetting either 8 or 12 identical volumes into titration plates
- ◆ Designed for operator comfort over long term repetitive pipetting
- ◆ Viton ring sealing provides tight fitting of tips
- ◆ One-handed adjustment and 4-digit display with prismatic lens window for easy viewing

### Transferpette® S - 8

8-channel pipettor. Supplied complete with performance certificate, a pack of pipette tips in a rack, a shelf/wall rack hanger, a set of Viton shaft seals and a reagent reservoir.

	Range µl	Subdiv. µl	A%±	V%	Tips reqd.
<b>PR174-35</b>	0.5 to 10	0.05	1.6	1.0	Crystal
<b>PR174-45</b>	5 to 50	0.1	0.8	0.4	Yellow
<b>PR174-55</b>	10 to 100	0.1	0.8	0.3	Yellow/Natural
<b>PR174-60</b>	20 to 200	1.0	0.8	0.3	Yellow/Natural
<b>PR174-65</b>	30 to 300	1.0	0.6	0.3	Yellow

### Transferpette® S - 12

12-channel pipettor. Supplied complete with performance certificate, a pack of pipette tips in a rack, a shelf/wall rack hanger, a set of Viton shaft seals and reagent reservoir.

	Range µl	Subdiv. µl	A%±	V%	Tips reqd.
<b>PR174-75</b>	0.5 to 10	0.05	1.6	1.0	Crystal
<b>PR174-80</b>	5 to 50	0.1	0.8	0.4	Yellow
<b>PR174-85</b>	10 to 100	0.1	0.8	0.3	Yellow/Natural
<b>PR174-90</b>	20 to 200	1.0	0.8	0.3	Yellow/Natural
<b>PR174-95</b>	30 to 300	1.0	0.6	0.3	Yellow

A% = Accuracy, V% = Reproducibility.

Figures quoted are based on the maximum volume and distilled water at 20%.

- ◆ Central pipetting and independent tip ejector buttons are ergonomically positioned for operation in either hand and when wearing gloves
- ◆ All controls clearly separated to avoid errors
- ◆ Fully autoclavable (at 121° for 20 minutes)
- ◆ Manifold can be rotated for ease of use when pipetting a titration plate
- ◆ Easy calibration without tools

### Disposable Pipette Tips, Brand

High quality polypropylene tips which are Conformity Certified when used with Transferpette® pipettors. Cadmium free. Also suitable for use with Eppendorf and Gilson pipettors.

### TipBox Set

Supplied as a steam sterilisable plastic box with hinged transparent cover containing a filled tip rack, plus additional tip racks as indicated.

### TipRacks

Refill packs for TipBox Sets, comprising racks for tips to be dropped into empty TipBoxes as required.

### Loose Packs

Supplied loose in a bag or carton in the quantities indicated.

### Crystal, for volumes 0.5 to 20µl

**PR140-12** TipBox set, 5 racks x 96 tips

**PR140-14** TipRacks for PR140-12, 10 racks x 96 tips

**PR140-16** Loose pack of 2000 tips

### Yellow, for volumes 2 to 200µl

**PR140-24** TipBox set, 5 racks x 96 tips

**PR140-26** TipRacks for PR140-24, 10 racks x 96 tips

**PR140-28** Loose pack of 1000 tips

### Natural\*, for volumes 50 to 1000µl

**PR140-45** TipBox set, 5 racks x 96 tips

**PR140-48** TipRacks for PR140-45, 10 racks x 96 tips

**PR140-52** Loose pack of 1000 tips

\*Except for PR140-52 which are blue coloured tips.

### Natural, for volumes 500 to 5000µl

**PR140-80** Loose pack of 200 tips

### Natural, for volumes 1000 to 10000µl

**PR140-95** Loose pack of 200 tips



PR174-35 to PR174-65



PR174-75 to PR174-95

### Conformity Certified

Indicated by the symbol  $\text{H}$  shows that the product has been manufactured by Brand under quality control procedures according to Eichordnung the German Federal Weights and Measures Regulations.

Provided the instrument is used with the recommended accessories and in accordance with the supplied instructions, Brand certify conformity with annex 12 of the Eichordnung.