Tissue Culture

Freezer Racks

Chest Freezer Racks

Stainless steel. With handle and retaining bar (to secure boxes during transportation). Dimensions stated are overall W x D x H. Supplied singly and empty.

	Dimensions, mm	Shelves	
TP466-32	150 x 160 x 539	8	
TP466-36	150 x 160 x 605	9	

Stainless steel. For use particularly with cryogenic boxes up to maximum dimensions 135 x 135 x 51mm L x W x H. With handle and retaining bar (to secure boxes during transportation). Dimensions stated are overall W x D x H. Supplied singly and empty.

	Dimensions, mm	Shelves	
RG416-90	140 x 140 x 335	6	
RG416-95	140 x 140 x 610	11	

Note

These racks are not freezer-specific and therefore we suggest that dimensions are checked before ordering.

Inventory Systems, Lab Tower

Flexible, multi-deck storage rack/inventory tower systems to hold 1.5ml conical micro tubes, cryogenic vials and other tubes up to 12mm o.d. with maximum height 55mm, in removable racks of 81 (9 x 9 format). Ideal for benchtop transport, or freezer storage, and wet/dry incubation. Polycarbonate construction with very wide operating temperature range, -196° C to $+121^{\circ}$ C. The small 16 x 16cm footprint minimises valuable space requirements. The internal dimensions of each deck are 142 x 142mm and the standard rack can be replaced by any storage rack, tray or box up to these dimensions and maximum 60mm high. Each deck platform has a front lip to secure the rack or box and the top deck has a carry-grip* mounted on it to aid portability. Available as standard tower systems comprising 5, 8 or 10 decks with or without standard racks, or alternatively as components for self assembly by the user.

Inventory Systems, Lab Tower, assembled

As described. The dimension stated is overall assembled height. All systems are 160 x 160mm W x D. Supplied singly.

	Decks	With racks	Height, mm	
TP468-05	5	\checkmark	340	
TP468-07	5	Х	340	
TP468-16	8*	1	660	
TP468-18	8*	х	660	
TP468-27	10*	1	820	
TP468-29	10*	х	820	

*The 8 and 10-deck systems are supplied with 2 carry-grips to allow "splitting" of the towers.

Inventory System components

TP468-72	Lab Tower decks (supplied with two carry-grips), pack of
TP468-76	81-place polycarbonate tube racks, pack of 10
TP468-78	Carry-grips, pack of 5

Dry-ice makers - see IC410.


Igloo ice containers and trays - see IC475.

Storage boxes - see TP460/TP462.

Ultralow freezers - see RG400/RG415.



TP468 in use



TP466, RG416 are similar



© 0808 200 5227 e-mail: sales@thelabwarehouse.com fax orders: 020 7729 2657



16